

## Product Datasheet

# Mouse Prolactin Receptor Recombinant Protein (orb1819912)

<b>Description</b>	The Mouse Prolactin Receptor yeast-derived recombinant protein is not tagged and is naturally endoto...
<b>Target</b>	Prolactin Receptor
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QSPPGKPEIH KCRSPDKETF TCWWNPGSDG GLPTNYSLTY SKEGEKNTYE CPDYKTS GPN SCFFSKQYTS IWKIYIITVN ATNEMGSSTS DPLYVDVTYI VEPEPPRNLTV LEVKQLKDKK TYLWVKWLPP TITDVKTGWF TMEYEIRLKS EEADWEIHF TGHQTQFKVF DLYPGQKYL V QTRCKPDHGY WSRWGQEKSI EIPNDFTLKD (210)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	24.6 kDa
<b>Entrez</b>	<a href="#">19116</a>
<b>Expiration Date</b>	6 months from date of receipt.