

Product Datasheet

Mouse Prolactin Receptor Recombinant Protein (orb1819914)

Description	The Mouse Prolactin Receptor yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	Prolactin Receptor
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QSPPGKPEIH KCRSPDKETF TCWWNPGSDG GLPTNYSLTY SKEGEKNTYE CPDYKTSGPN SCFFSKQYTS IWKIYIITVN ATNEMGSSTS DPLYVDVTYI VEPEPPRNLTV LEVKQLKDKK TYLWVKWLPP TITDVKTGWF TMEYEIRLKS EEADWEIHF TGHQTQFKVF DLYPGQKYLTV QTRCKPDHGY WSRWGQEKSI EIPNDFTLKD (210)
Purity	98%
Source	Yeast
MW	24.6 kDa
Entrez	19116
Expiration Date	6 months from date of receipt.