

Product Datasheet

Human Prolactin Recombinant Protein (orb1819916)

Description	The Human Prolactin yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	Prolactin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LPICPGGAAR CQVTLRDLFD RAVVLSHYIH NLSSEMFSEF DKRYTHGRGF ITKAINSCHT SSLATPEDKE QAQQMNQKDF LSLIVSILRS WNEPLYHLVT EVRGMQEAPE AILSKAVEIE EQTKRLLEGM ELIVSQVHPE TKENEIYPVW SGLPSLQ MAD EESRLSAYYN LLHCLRRDSH KIDNYLKKL CRIIHNNNC (199)
Purity	98%
Source	Yeast
MW	22.9 kDa
Entrez	5617
Expiration Date	6 months from date of receipt.