

Product Datasheet

Human Osteopontin Recombinant Protein (orb1819922)

Description	The Human Osteopontin yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	Osteopontin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IPVKQADSGS SEEKQLYNKY PDAVATWLNPDPSQKQNLLA PQNAVSSEET NDFKQETLPS KSNESHDMMDMDDEDDDDH VDSQDSIDSNDSDDVDDTDD SHQSDSHHSDESDELVTDF PTDLPATEVFTPVVPTVDY DGRGDSVYGLRSKSKKFRF PDIQYPDATDEDITSHMESE ELNGAYKAIPVAQDLNAPSD WDSRGKDSYETSQDDQSAE THSHKQSRLYKRKANDESNE HSDVIDSQELSKVSREFHSH EFHSHEDMLV VDPKSKEEDK HLKFRISHELDSASSEVN (298)
Purity	98%
Source	Yeast
MW	33.7 kDa
Entrez	6696
Expiration Date	6 months from date of receipt.