

Product Datasheet

Bovine PD-L1 Recombinant Protein (orb1819941)

Description	The Bovine PD-L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	PD-L1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FTITVSKDLY VVEYGSNVTLECRFPVDKQL NLLVLVYWE MEDKKIIQFV NGKEDPNVQH SSYHGRAQLL KDQLFLGKAA LQITDVKLQD AGVYCLISY GGADYKRITL KVNAPYRKIY HTISVDPVTS EHELTCQAEG YPEADVIWTS SDHQVLSGKT SITSSKREEK LFNVTSTLRI NTTADKIFYC TFRRLGHEEN NTAELVIPEP YLDPAKRRNH (220)
Purity	98%
Source	Yeast
MW	25.1 kDa
Entrez	533834
Expiration Date	6 months from date of receipt.