

Product Datasheet

Human CCL20 Recombinant Protein (orb1819952)

Description	The Human CCL20 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL20
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ASNFDCCCLGY TDRILHPKFI VGFTRQLANE GCDINAIIFH TKKKLSVCAN PKQTWVKYIV RLLSKKVKNM (70)
Purity	98%
Source	Yeast
MW	8.0 kDa
Entrez	6364
Expiration Date	6 months from date of receipt.