

Product Datasheet

Turkey CCL3 Recombinant Protein (orb1819953)

Description	The Turkey CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IKVAAQRPS IVQSRALQPP SSSPAHCANV LQTMKGPAAI LTALLLLALC SSAVAHLDEL PTTCCFSYVQ RPIPRNIIKS AYITSSKRL PAVILVTKKG KEICANPEES WVQNHLELFQ NEEN (124)
Purity	98%
Source	Yeast
MW	13.5 kDa
Entrez	100544108
Expiration Date	6 months from date of receipt.