

Product Datasheet

Caprine CCL3 Recombinant Protein (orb1819956)

Description	The Caprine CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APFGADTPTA CCFSYAARQL SRKIIADYFE TSSQCSKPGV IFQTKRGRQV CANPSEDWVQ EYVTDLELNS (70)
Purity	98%
Source	Yeast
MW	7.8 kDa
Entrez	102181154
Expiration Date	6 months from date of receipt.