

## Product Datasheet

### Human CCL18 Recombinant Protein (orb1819965)

|                         |   |
|-------------------------|---|
| <b>Description</b>      | The Human CCL18 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| <b>Target</b>           | CCL18   |
| <b>Form/Appearance</b>  | Lyophilized   |
| <b>Storage</b>          | -20°C   |
| <b>Note</b>             | For research use only   |
| <b>Protein Sequence</b> | AQVGTNKELC CLVYTSWQIP QKFIVDYSET<br>SPQCPKPGVI LLTKRGRQIC<br>ADPNKKWVQK YISDLKLNA (69)                  |
| <b>Purity</b>           | 98%   |
| <b>Source</b>           | Yeast   |
| <b>MW</b>               | 7.9 kDa   |
| <b>Entrez</b>           | <a href="#">6362</a>  |
| <b>Expiration Date</b>  | 6 months from date of receipt.  |