

Product Datasheet

Human CCL18 Recombinant Protein (orb1819966)

Description	The Human CCL18 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL18
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AQVGTNKELC CLVYTSWQIP QKFIVDYSET SPQCPKPGVI LLTKRGRQIC ADPNKKWVQK YISDLKLNA (69)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	6362
Expiration Date	6 months from date of receipt.