

Product Datasheet

Human Amphiregulin Recombinant Protein (orb1819970)

Description	The Human Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	Amphiregulin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SVRVEQVVKP PQNKTESENT SDKPKRKKKG GKNGKNRRNR KKKNPCNAEF QNFCIHGECK YIEHLEAVTC KCQQEYFGER CGEKSMKTHS MIDSSLSK (98)
Purity	98%
Source	Yeast
MW	11.3 kDa
Entrez	374
Expiration Date	6 months from date of receipt.