

Product Datasheet

Human CCL3 Recombinant Protein (orb1819976)

| | |
|-------------------------|---|
| Description | The Human CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ... |
| Target | CCL3 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | ASLAADTPTA CCFSYTSRQI PQNFIADYFE TSSQCSKPGV IFLTKRSRQV CADPSEEWVQ KYVSDLELSA (70) |
| Purity | 98% |
| Source | Yeast |
| MW | 7.8 kDa |
| Entrez | 6348 |
| Expiration Date | 6 months from date of receipt. |