

Product Datasheet

Mouse TIM-3 Recombinant Protein (orb1819982)

Description	The Mouse TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	TIM-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RSLENAYVFE VGKNAYLPCS YTLSTPGALV PMCWKGKFCP WSQCTNELLR TDERNVTYQK SSRYQLKGD L NKGDVSLIK NVTLDDHGT Y CCRIQFPGLM NDKKLELKLD IKA AKVTPAQ TAHGDSTTAS PRTL TTERN G SETQTLVTLH N NNGTKISTW ADEIKDSGET IRTA (174)
Purity	98%
Source	Yeast
MW	19.3 kDa
Entrez	171285
Expiration Date	6 months from date of receipt.