

Product Datasheet

Mouse TIM-3 Recombinant Protein (orb1819983)

| | |
|-------------------------|--|
| Description | The Mouse TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | TIM-3 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | RSLENAYVFE VGKNAYLPCS YTLSTPGALV PMCWGKGFPC WSQCTNELLR TDERNVTYQK SSRYQLKGD L NKGDVSLIK NVTLDDHGT Y CCRIQFPGLM NDKKLELKLD IKA AKVTPAQ TAHGDSTTAS PRTL TTERN G SETQTLVTLH N NNGTKISTW ADEIKDSGET IRTA (174) |
| Purity | 98% |
| Source | Yeast |
| MW | 19.3 kDa |
| Entrez | 171285 |
| Expiration Date | 6 months from date of receipt. |