

Product Datasheet

Bovine PlGF Recombinant Protein (orb1819987)

Description	The Bovine Placental Growth Factor yeast-derived recombinant protein is not tagged and is naturally ...
Target	PlGF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LPAVPTTQWA LSPGNISSEV EVVPFQQVWS RSYCRPVERL VDIVSEYPSE MEHLFSPSCV SLMRCTGCCS DESMHCVPLE TANVTMQLMK YRSLDQPPFV EMSFSQHVRC ECKPLWEKMK QTRCGDTISQ R (131)
Purity	98%
Source	Yeast
MW	15.1 kDa
Entrez	280894
Expiration Date	6 months from date of receipt.