

Product Datasheet

Rabbit IL-16 Recombinant Protein (orb1819994)

Description	The Rabbit IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GLSPSSDSAA APSAASDVPA ESATEAAICT VTLEKTSAGL GFSLEGGRGS LRGDKPITIN RIFKGVASGP SATLQPGDEI LHVAGTALQG LTRFEAWNII KALPDGPVTV VIRRLQPQG TVATGDP (127)
Purity	98%
Source	Yeast
MW	12.7 kDa
Entrez	100356333
Expiration Date	6 months from date of receipt.