

Product Datasheet

Atlantic Salmon IL-8 Recombinant Protein (orb1819998)

Description	The Atlantic Salmon IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MSLRGVGADL RCRCIETESR RIGRLIKKVE MFPPSSHCRD TEIIATLSNS GQEICLDVSA PWVKRVIKKM LANNK (75)
Purity	98%
Source	Yeast
MW	8.5 kDa
Entrez	100195681
Expiration Date	6 months from date of receipt.