

## Product Datasheet

### Chicken EGF Recombinant Protein (orb1820009)

<b>Description</b>	The Chicken EGF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	EGF
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Application notes</b>	Chicken EGF applications are for cell culture, ELISA standard, and Western Blot Control.
<b>Protein Sequence</b>	GDSIGCPPAY DSYCLHGGVC NYVSDLQDYA CNCVTGYVGE RCQFSDLEWW EQQHA (55)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	6.2 kDa
<b>Entrez</b>	<b>408035</b>
<b>Expiration Date</b>	6 months from date of receipt.