

Product Datasheet

Rabbit IFN-beta Recombinant Protein (orb1820027)

Description	The Rabbit IFN beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IFN beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MSYNSLQIQL WHGSLTCAKL LLQLNGTTED CLNERINFKV PKEIKEPQQL QKEDTTLVIF EMLNNIFDIF RKNFSSTGWN ETLVENLLGE THLQIHHLKS KINKKVTLES IRMNLRKSY YWRIMDYLET KQYSNCAWKI VQLEIFRNFS FIIMLIDYL (159)
Purity	98%
Source	Yeast
MW	19.0 kDa
Entrez	100342931
Expiration Date	6 months from date of receipt.