

Product Datasheet

Chicken LIF Recombinant Protein (orb1820028)

Description	The Chicken Leukemia Inhibitory Factor yeast-derived recombinant protein is not tagged and is natura...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RALLGTSSAC PTNGLCRANV LEQTRRQVAL LNATAQDLFS LYLKCQGEPE SSESRLCSP SGIFFPPFHV NRTTERKEVM VAMYKLF AFL NASLGNITRD QEELNPMKAE LLNRLHNTTK TTRGLISNLT CLLCKHYNIF QVDVSYGESS KDKSAFKKKQ QGCQVLRKYV QVIAQAARVL LPHLSPA (187)
Purity	98%
Source	Yeast
MW	21.0 kDa
Entrez	427718
Expiration Date	6 months from date of receipt.