

## Product Datasheet

# Mouse PDGF-BB Recombinant Protein (orb1820035)

<b>Description</b>	The Mouse PDGF-BB yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
<b>Target</b>	PDGF-BB
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SLGSLAAAEP AVIAECKTRT EVFQISRNLI DRTNANFLVW PPCVEVQRCS GCCNNRNVQC RASQVQMRPV QVRKIEIVRK KPIFKKATVT LEDHLACKCE TIVTPRPVT (109)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	12.2 kDa (monomer)
<b>Entrez</b>	<a href="#">18591</a>
<b>Expiration Date</b>	6 months from date of receipt.