

Product Datasheet

Mouse PDGF-BB Recombinant Protein (orb1820037)

Description	The Mouse PDGF-BB yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	PDGF-BB
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SLGSLAAAEP AVIAECKTRT EVFQISRNLI DRTNANFLVW PPCVEVQRCS GCCNNRNVQC RASQVQMRPV QVRKIEIVRK KPIFKKATVT LEDHLACKCE TIVTPRPVT (109)
Purity	98%
Source	Yeast
MW	12.2 kDa (monomer)
Entrez	18591
Expiration Date	6 months from date of receipt.