

Product Datasheet

Human IL-18 Recombinant Protein (orb1820039)

Description	The Human IL-18 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-18
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	YFGKLESKLS VIRNLNDQVL FIDQGNRPLF EDMTDSDCRD NAPRTIFIIS MYKDSQPRGM AVTISVKCEK ISTLSCENKI ISFKEMNPPD NIKDTKSDII FFQRSVPGHD NKMQFESSY EGYFLACEKE RDLFKLILKK EDELGDRSIM FTVQNE (157)
Purity	98%
Source	Yeast
MW	18.2 kDa
Entrez	3606
Expiration Date	6 months from date of receipt.