

Product Datasheet

Atlantic Salmon CCL4 Recombinant Protein (orb1820050)

Description	The Atlantic Salmon CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-...
Target	CCL4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LRMASGPKKC CTFADRQIP RGKVVGYYTKT SQQCSNPAIL FKTQAGRQVC ARSSDRWVKD YINFLDSKKS GEQTPLL (77)
Purity	98%
Source	Yeast
MW	8.7 kDa
Entrez	100136511
Expiration Date	6 months from date of receipt.