

Product Datasheet

Atlantic Salmon CCL4 Recombinant Protein (orb1820051)

| | |
|-------------------------|--|
| Description | The Atlantic Salmon CCL4 yeast-derived recombinant protein is not tagged and is naturally endotoxin... |
| Target | CCL4 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | LRMASGPKKC CTFADRQIP RGKVVGYTKT SQQCSNPAIL FKTQAGRQVC ARSSDRWVKD YINFLDSKKS GEQTPLL (77) |
| Purity | 98% |
| Source | Yeast |
| MW | 8.7 kDa |
| Entrez | 100136511 |
| Expiration Date | 6 months from date of receipt. |