

Product Datasheet

Atlantic Salmon IL-6 Recombinant Protein (orb1820055)

Description	The Atlantic Salmon IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPVPSALSEL MTSGWTSGEE LATDSLYGET GAPPKWEQII KMLVHEVTTL RNQQFVEEFQ KPVVEISSFS QHQVLSTPPH LSKTLCSTTN KEACLQEIRR GLRVYQLLLQ HVKAEYPQST LLPSVTHQTT VLIQLVKDQM KAAEVVEDLS ASERERVLGE VSTGTEWERK TSVHAILREL RNFLVDTKRA LRRMEKRGKD FQ (202)
Purity	98%
Source	Yeast
MW	22.9 kDa
Entrez	106570581
Expiration Date	6 months from date of receipt.