

Product Datasheet

Human IL-13 Recombinant Protein (orb1820056)

Description	The Human IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPVPPSTALRELIEELVNITQNQKAPLCNGSMVWSINLTAGMYCAALESLINVSG CSAIEKTQRMLSGFCPHKVSAGQFSSLHVRDTKIEVAQFVKDLLLHLKLFREGQ FN (112)
Purity	98%
Source	Yeast
MW	12.3 kDa
Entrez	3596
Expiration Date	6 months from date of receipt.