

Product Datasheet

Human IL-13 Recombinant Protein (orb1820057)

| | |
|-------------------------|---|
| Description | The Human IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The... |
| Target | IL-13 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | GPVPPSTALRELIEELVNITQNQKAPLCNGSMVWSINLTAGMYCAALESLINVSG CSAIEKTQRMLSGFCPHKVSAGQFSSLHVRDTKIEVAQFVKDLLLHLKLFREGQ FN (112) |
| Purity | 98% |
| Source | Yeast |
| MW | 12.3 kDa |
| Entrez | 3596 |
| Expiration Date | 6 months from date of receipt. |