

## Product Datasheet

### IGF-1 Recombinant Protein (orb1820064)

<b>Description</b>	The IGF-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The purit...
<b>Target</b>	IGF-1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GPETLCGAEL VDALQFVCGD RGFYFNKPTG YGSSRRAPQ TGIVDECCFR SCDLRRLEMY CAPLKPAKSA (70)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	7.6 kDa
<b>Expiration Date</b>	6 months from date of receipt.