

Product Datasheet

IGF-1 Recombinant Protein (orb1820065)

Description	The IGF-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The purit...
Target	IGF-1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPETLCGAEL VDALQFVCGD RGFYFNKPTG YGSSRRAPQ TGIVDECCFR SCDLRRLEMY CAPLKPAKSA (70)
Purity	98%
Source	Yeast
MW	7.6 kDa
Expiration Date	6 months from date of receipt.