

## Product Datasheet

# Atlantic Salmon TNF alpha Recombinant Protein (orb1820073)

<b>Description</b>	The Atlantic Salmon TNF alpha yeast-derived recombinant protein is not tagged and is naturally endot...
<b>Target</b>	TNF alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AAIHLEGEYN PNLTADTVQW RKDDGQAFSQ GGFELQGNQI LIPHTGLFFV YSQASFRVKC NGPGENTTPL SHVIWRYSDS IGVNANLLSG VRSVCQRNYG DAESNIGEAW YNAVYLGAVF QLNEGDKLWT ETNRLTNVEP EQGKNFFGVF EL (152)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.9 kDa
<b>Entrez</b>	<a href="#">100136457</a>
<b>Expiration Date</b>	6 months from date of receipt.