

Product Datasheet

Atlantic Salmon TNF alpha Recombinant Protein (orb1820074)

Description	The Atlantic Salmon TNF alpha yeast-derived recombinant protein is not tagged and is naturally endot...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAIHLEGEYN PNLTADTVQW RKDDGQAFSQ GGFELQGNQI LIPHTGLFFV YSQASFRVKC NGPGENTTPL SHVIWRYSDS IGVNANLLSG VRSVCQRNYG DAESNIGEAW YNAVYLGAVF QLNEGDKLWT ETNRLTNVEP EQGKNFFGVF EL (152)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	100136457
Expiration Date	6 months from date of receipt.