

Product Datasheet

Canine PDGF-BB Recombinant Protein (orb1820085)

Description	The Canine PDGF-BB yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	PDGF-BB
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SLGSPTVAEP AMIAECKTRT EVFEISRRLI DRTNANFLVW PPCVEVQRCS GCCNNRNVQC RPTQVQLRPV QVRKIEIVRK KPTFKKATVT LEDHLACKCE TVVAARPVT (109)
Purity	98%
Source	Yeast
MW	12.3 kDa (monomer)
Entrez	442986
Expiration Date	6 months from date of receipt.