

Product Datasheet

Zebrafish TNF alpha Recombinant Protein (orb1820091)

Description	The Zebrafish TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
Target	TNF alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AAIHLTGGYN SESKTLDWRD DQDQAFSSGG LKLVNREIII PDDGIYFVYS QVSLHISCTS ELTEEQVLMS HAVMRFSESY GGKKPLFSAI RSICTQEPES ENLWYNTIYL GAAFHLREGD RLGTDTTTAL LPMVENDNGK TFFGVFGL (148)
Purity	98%
Source	Yeast
MW	16.5 kDa
Entrez	405785
Expiration Date	6 months from date of receipt.