

## Product Datasheet

### Cannine PD-L1 Recombinant Protein (orb1820095)

<b>Description</b>	The Canine PD-L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	PD-L1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	FTITVSKDLY VVEYGGNVTM ECKFPVEKQL NLFALIVYWE MEDKKIIQFV NGKEDLKVQH SSYSQRAQLL KDQLFLGKAA LQITDVRLQD AGVYCCLIGY GGADYKRITL KVHAPYRNIS QRISVDPVTS EHELMCQAEG YPEAEVIWTS SDHRVLSGKT TITNSNREEK LFNVTSTLNI NATANEIFYC TFQRSGPEEN NTAELVIPER LPVPASER (218)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	24.7 kDa
<b>Entrez</b>	<a href="#">484186</a>
<b>Expiration Date</b>	6 months from date of receipt.