

Product Datasheet

Cannine PD-L1 Recombinant Protein (orb1820097)

Description	The Canine PD-L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	PD-L1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FTITVSKDLY VVEYGGNVTM ECKFPVEKQL NLFALIVYWE MEDKKIIQFV NGKEDLKVQH SSYSQRAQLL KDQLFLGKAA LQITDVRLQD AGVYCCLIGY GGADYKRITL KVHAPYRNIS QRISVDPVTS EHELMCQAEG YPEAEVIWTS SDHRVLSGKT TITNSNREEK LFNVTSTLNI NATANEIFYC TFQRSGPEEN NTAELVIPER LPVPASER (218)
Purity	98%
Source	Yeast
MW	24.7 kDa
Entrez	484186
Expiration Date	6 months from date of receipt.