

Product Datasheet

Mouse SCF Recombinant Protein (orb1820103)

| | |
|-------------------------|--|
| Description | The Mouse SCF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p... |
| Target | SCF |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | KEICGNPVTD NVKDITKLVA NLPNDYMITL NYVAGMDVLP SHCWLRDMVI QLSLSLTLL DKFSNISEGL SNYSIIDKLG KIVDDLVLCM EENAPKNIKE SPKRPETRSF TPEEFSIFN RSIDAFKDFM VASDTSDCVL SSTLGPEKDS RVSVTKPFML PPVAA (165) |
| Purity | 98% |
| Source | Yeast |
| MW | 18.4 kDa |
| Entrez | 17311 |
| Expiration Date | 6 months from date of receipt. |