

Product Datasheet

Swine IL-3 Recombinant Protein (orb1820107)

| | |
|-------------------------|--|
| Description | The Swine IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ... |
| Target | IL-3 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | MPTTTLQPKN YLAMIQEITR SLENLTVTSN KSLTLNELET LVNNTLLRPN LEAFVTF AEN HLKNISGIKK NLEKFRPILP TSMSTEEPIS IEEGDLGDFR AKLMEYLVVL RDSLKPMITE P (121) |
| Purity | 98% |
| Source | Yeast |
| MW | 13.7 kDa |
| Expiration Date | 6 months from date of receipt. |