

Product Datasheet

Rabbit IL-6 Recombinant Protein (orb1820121)

Description	The Rabbit IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPTSAPVRED SNTKASPKT LTPPGRTIES IRSILETIKE LRKEMCDHDV NCMNRKEALA EVNLHLPRLI EEDGCFPPAV NNETCLLRIT SGLMEFRMYL EHLQAKFRSD EENTRVSMVL KNIQHLLIKTL RPKVKNLNEE ATLKPAVAVS LMENLQQKNQ WLKTTTIHFI LRGLTNFLEF TLRAVDLM (188)
Purity	98%
Source	Yeast
MW	24.5 kDa
Entrez	100008733
Expiration Date	6 months from date of receipt.