

Product Datasheet

IRF-4 Recombinant Protein (orb1820122)

Description	The IRF-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The purit...
Target	IRF-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NGKLRQWLID QIDSGKYPGL VWENEEKSIF RIPWKHAGKQ DYNREEDAAL FKAWALFKGK FREGIDKPDP PTWKTRLRCA LNKSNDFEEL VERSQLDISD PYKVYRIVP (109)
Purity	98%
Source	Yeast
MW	12.9 kDa
Expiration Date	6 months from date of receipt.