

Product Datasheet

FNDC-5/Irisin Recombinant Protein (orb1820126)

Description	The FNDC5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The purit...
Target	FNDC-5
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DSPSAPVNVV VRHLKANSV VSWDVLEDEV VIGFAISQK KDVRMLRFIQ EVNTTTRSCA LWDLEEDTEY IVHVQAISIQ GQSPASEPVL FKTPREAEM ASKNKDEVTM KE (112)
Purity	98%
Source	Yeast
MW	12.6 kDa
Expiration Date	6 months from date of receipt.