

Product Datasheet

FNDC-5/Irisin Recombinant Protein (orb1820127)

Description The FNDC5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The purit...

Target FNDC-5

Form/Appearance Lyophilized

Storage -20°C

Note For research use only

Protein Sequence DSPSAPVNVV VRHLKANSV
VSWDVLEDEV VIGFAISQK
KDVRMLRFIQ EVNTTTRSCA
LWDLEEDTEY IVHVQAISIQ GQSPASEPVL
FKTPREAEM ASKNKDEVTM KE (112)

Purity 98%

Source Yeast

MW 12.6 kDa

Expiration Date 6 months from date of receipt.