

Product Datasheet

Bovine IL-3 Recombinant Protein (orb1820128)

Description	The Bovine IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APQAKGLPVV TS RTPYSMLM KEIMDDLKKI TPSPEGLNS DEKNFLTKES LLQANLKVFM TFATDTFGSD SKIMKNLKEF QPVLPTATPT EDPIFIENKN LGDFRMKLEE YLVIIRNYLK SKNIWFS (127)
Purity	98%
Source	Yeast
MW	14.5 kDa
Entrez	280823
Expiration Date	6 months from date of receipt.