

Product Datasheet

Rabbit PD-L1 Recombinant Protein (orb1820133)

Description	The Rabbit PD-L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	PD-L1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FTVTVPKDLY VVEYGSNVTI ECKFPVEKQL DLNSLVVYWE MGDKQIIQFV HG NENLKVQH SSYRQRAQLL KDQLHLGIAA LQITDVKLRD AGVYCCMISY GGADYKRITL KVNAPYHKIN QRISVDPVTS EHELTCQAEG YPEAEVIWTS SDHQILHGKT DVTNSNREEK FFNVTSTLRI NTTANEIFYC TFQRPGPEAN NTAKLIIPEP PVFHRNER (218)
Purity	98%
Source	Yeast
MW	24.9 kDa
Entrez	100357106
Expiration Date	6 months from date of receipt.