

Product Datasheet

Canine EGF Recombinant Protein (orb1820137)

Description	The Canine EGF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	EGF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NGYRECPSSY DGYCLYNGVC MYIEAVDRYA CNCVFGYVGE RCQHRDLKWE LR (52)
Purity	98%
Source	Yeast
MW	6.2 kDa
Entrez	403657
Expiration Date	6 months from date of receipt.