

Product Datasheet

Rabbit IL-1 Receptor Antagonist Recombinant Protein (orb1820141)

Description	The Rabbit IL-1 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturally...
Target	IL-1RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SGKRPCRMQA FRIWDVNQKT FYLRNNQLVA GYLQGPNAKL EERIDVVPLE PQLFLGIQR GKLCLSCVKS GDKMKLHLEA VNITDLGKNK EQDKRFTFIR SNSGPTTTFE SASCPGWFLC TALEADQPVS LTNTPDDSI VTKFYFQEDQ (150)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	100009558
Expiration Date	6 months from date of receipt.