

Product Datasheet

Rabbit IL-13 Recombinant Protein (orb1820174)

Description	The Rabbit IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPVPPPTSLK ELIEELVNIT HNQKAPLCNG TMVWSVNLTG SVYCAALES VNVSGCNAIQ RTQRMLSGLC TDKAVAKQVT SVQARDTKIE LLQFLKELRR HLQMLYRLGK FR (112)
Purity	98%
Source	Yeast
MW	12.5 kDa
Entrez	100358676
Expiration Date	6 months from date of receipt.