

Product Datasheet

Rat IL-4 Recombinant Protein (orb1820186)

Description	The Rat IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The pu...
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CNDSPLREII NTLNQVTEKG TPCTEMFVPD VLTATRNTTE NELICRASRV LRKFYFPRDV PPCLKNKSGV LGELRKLCRG VSGLNSLRSC TVNESTLTTL KDFLESLSKI LRGKYLQSCT SMS (123)
Purity	98%
Source	Yeast
MW	13.8 kDa
Entrez	287287
Expiration Date	6 months from date of receipt.