

Product Datasheet

Dolphin CXCL10 Recombinant Protein (orb1820193)

Description	The Dolphin CXCL10 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	CXCL10
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IPLSRTTRCT CIKISDRPVN PRSLEKLEVI PASQSCPRVE IIATMKKNGE KRCLNPESKT IKNLLKAISK ERSKRSQTQK EA (82)
Purity	98%
Source	Yeast
MW	9.3 kDa
Entrez	101336308
Expiration Date	6 months from date of receipt.