

Product Datasheet

Ferret IFN gamma Recombinant Protein (orb1820196)

Description	The Ferret IFN gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CDLPQDHGLL HWRALMLLQQ MRRLSASSCD KYTNDFGFPQ EVFDGKQLQK AQALSVIHVT NQKTFHLFCT EASPAPWNTT LLEQLCSGLS EQLGHAAAC LQEAGVGELP LVNGDSILRN YFQRISLYLQ EKQYSPCAWE MVRAEIMKPL YASTVLHKRL RSRK (164)
Purity	98%
Source	Yeast
MW	18.7 kDa
Entrez	101693741
Expiration Date	6 months from date of receipt.